

# Recombinant Human Protein FAM172A,FAM172A (C-6His)



Catalog No: #AP76006

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP76006-1 10ug #AP76006-2 50ug #AP76006-3 500ug #AP76006-4 1mg

## Description

Product Name	Recombinant Human Protein FAM172A,FAM172A (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q8WUF8
Uniprot	Q8WUF8
GeneID	83989;
Calculated MW	47kDa
Target Sequence	QIQQGGPDEKEKTTALKDLLSRIDLDELMMKKDEPPLDFPDTLEGFEYAFNEKGQLRHIKTGEPFVFNYREDLHR WNQKRYEALGEIITKYVYELLEKDCNLKKSIPVDATSEPKSFIFMSEDALTNPQKLMVLIHGSGVVRAGQWA RRLIINEDLDSGTQIPFIKRAVAEGYGVIVLNPENYIEVEKPKIHVQSSSDSSDEPAEKREKDKVSKETKRR DFYEKYRNPQREKEMMQLYIRENGSPEEHAIYVWDHFIAQAAAENVFFVAHSYGGLAFVELMIQREADVKNKV TAVALTDSVHNVWHQEAGKTIREWMRENCCNWWSSSEPLDTSVESMLPDCPRVSAGTDRHELTSWKSFPSI FKFFTEASEAKTSSLKPAVTRRSHRIKHEELVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only