

# Recombinant Human Prostate-Specific Antigen,,PSA,KLK3 (Ile25-Pro261, C-6His)



Catalog No: #AP76009

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP76009-1 10ug #AP76009-2 50ug #AP76009-3 500ug #AP76009-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Prostate-Specific Antigen,,PSA,KLK3 (Ile25-Pro261, C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P07288
Uniprot	P07288
GeneID	354;
Calculated MW	27.1kDa
Target Sequence	IVGGWECEKHSQPWQVLVASRGRAVCGGVLVHPQWVLTAAHCIRNKSVILLGRHSLFHPEDTGQVFQVSHS FPHPLYDMSLLKNRFLRPGDDSSHDLMLLRLSEPAELTDAVKVMDLPTQEPALGTTTCYASGWGSIEPEEFLTP KKLQCVDLHVISNDVCAQVHPQKVKFMLCAGRWTGGKSTCSGDSGGPLVCNGVLQGITSWGSEPCALPER PSLYTKVVHYRKWIKDTIVANPVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHcl, 150mM NaCl, 1mM CaCl <sub>2</sub> , pH7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only