

Recombinant Human Cytotoxic T-Lymphocyte Protein 4,CTLA-4,CD152 (C-Fc)



Catalog No: #AP76010

Orders: order@signalwayantibody.com

Package Size: #AP76010-1 10ug #AP76010-2 50ug #AP76010-3 500ug #AP76010-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Cytotoxic T-Lymphocyte Protein 4,CTLA-4,CD152 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P16410
Uniprot	P16410
GeneID	1493;
Calculated MW	41.5kDa
Target Sequence	KAMHVAQPAVVLASSRGIASFVCEYASPGKATEVRVTVLRQADSQVTEVCAATYMMGNELTFLDDSICTGTSS GNQVNLTIQGLRAMDTGLYICKVELMYPPPYLIGINGTQIYVIDPEPCPDSQEPKSSDKTHTSPSPAPPELL GGSSVFLFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLT VLHQQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSRDELTKNQVSLTCLVKGFYPSDIAVE WESNGQPENNYKTTTPVLDSGSEFLYSLKLTVDKSRWQQGNVFCSSVMHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only