

Recombinant Human NPDC1,CAB1 (C-6His)

Catalog No: #AP76012



Package Size: #AP76012-1 10ug #AP76012-2 50ug #AP76012-3 500ug #AP76012-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human NPDC1,CAB1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9NQX5
Uniprot	Q9NQX5
GeneID	56654;
Calculated MW	16.5kDa
Target Sequence	GHPDVAACPGSLDCALKRRRARCPPGAHACGPCLQPFQEDQQGLCVPRMRRPPGGGRPQPRLEDEIDFLAQ ELARKESGHSTPPLPKDRQRLPEPATLGFSARGQGLELGLPSTPGTPTPTHTSLGSPVSSDPVHMSPLEPR GGQGDVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, pH8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only