

Recombinant Human Leucine-Rich Repeat-Containing Protein 25,LRRC25 (C-6His)

Catalog No: #AP76013

Orders: order@signalwayantibody.com

Package Size: #AP76013-1 10ug #AP76013-2 50ug #AP76013-3 500ug #AP76013-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Leucine-Rich Repeat-Containing Protein 25,LRRC25 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q8N386
Uniprot	Q8N386
GeneID	126364;
Calculated MW	16.7kDa
Target Sequence	LEPSTVSSADVDWNAEFSATCLNFGSLSLPHNQSLRASNVILLDLSGNGLRELPTVFFAHLQKLEVLNVLR NPLSRVDGALAARCDLDLQADCNCALESWHDIRRDNCSGQKPLLCWDTTSSQHNLSAFLEVSCAPGLASATV DHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only