

Recombinant Human α -Amylase 2B, AMY2B (C-6His)

Catalog No: #AP76014



Package Size: #AP76014-1 10ug #AP76014-2 50ug #AP76014-3 500ug #AP76014-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human α -Amylase 2B, AMY2B (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P19961
Uniprot	P19961
GeneID	280;
Calculated MW	56.9kDa
Target Sequence	QYSPNTQQGRTSIVHLFEWRWVDIALECERYLAPKGFGGVQVSPNENVAIHNPFRPWWERYQPVSYKLC RSGNEDEFNRNMVTRCANNVGVRIYVDAVINHMSGNAVSAGTSSTCGSYFNPGSRDFPAVPYSGWDFNDGKCK TGSGDIENYNDATQVRDCRLVGLLDLDALEKDYVRSKIAEYMNHLIDIGVAGFRLDASKHMWPGDIKAILDKLHN LNSNWFPAGSKPFIYQEVIDLGGEPIKSSDYFGNGRVTEFKYGAKLGTVIRKWNGEKMSYLKNWGEWGF PSDRALVFDNHDNQRGHGAGGASILTFWDARLYKMAVGFMLAHPYGFTRVMSSYRWPRQFQNGNDVND WVGPPNNGVIKEVTINPDTCGNDWVCEHRWRQIRNMVFRNVVDGQPFTNWYDNGSNQVAFGRGNRG FIVFNDDWTFSLTLQTGLPAGTYCDVISGDKINGNCTGIKIYVSDDGKAHFSISNSAEDPFIAIHAESKLVDDHH HHH
Formulation	Supplied as a 0.2 μ m filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only