

# Recombinant Mouse Placenta Growth Factor,PGF,PIGF (C-6His)



Catalog No: #AP76015

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP76015-1 10ug #AP76015-2 50ug #AP76015-3 500ug #AP76015-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                 |  |
|-----------------|--|
| Product Name    | Recombinant Mouse Placenta Growth Factor,PGF,PIGF (C-6His)   |
| Host Species    | Human Cells  |
| Purification    | Greater than 95% as determined by reducing SDS-PAGE.   |
| Specificity     | Mouse  |
| Accession No.   | P49764   |
| Uniprot         | P49764   |
| GeneID          | 18654;   |
| Calculated MW   | 16.9kDa  |
| Target Sequence | VHSQGALSAGNNSTEVEVVPFNEVWGRSYCRPMEKLVYILDEYPDEVSHIFSPSCVLLSRCSGCCGDEGLHC<br>VPIKTANITMQILKIPPNRDPHFYVEMTFSQDVLCECRPILETTKAERRKTKGKRKRSRNSQTEEPHPVDHHHH<br>HH   |
| Formulation     | Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.   |
| Storage         | Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months. |

Note: This product is for in vitro research use only