

Recombinant Mouse C-X3-C Motif Chemokine 1,CX3CL1,Fractalkine (C-6His)

Catalog No: #AP76016

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76016-1 10ug #AP76016-2 50ug #AP76016-3 500ug #AP76016-4 1mg

Description

Product Name	Recombinant Mouse C-X3-C Motif Chemokine 1,CX3CL1,Fractalkine (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	O35188
Uniprot	O35188
GeneID	20312;
Calculated MW	34.3kDa
Target Sequence	<p>QHLGMTKCEIMCDKMTSRIPVALLIRYQLNQESCGKRAIVLETTQHRRFCADPKEKWWQDAMKHLDHQAAALT</p> <p>KNGGKFEKRVNDNVTPGITLATRGLSPSALTKPESATLEDLALELTTISQEARGTMGTSQEPPAAVTGSSLSTSE</p> <p>AQDAGLTAKPQSIGSFEEADISTTVWPSPAVYQSGSSSWAEKATESPSTTAPSPQVSTTSPSTPEENVGSE</p> <p>GQPPWVQGQDLSPEKSLGSEEINPVHTDNFQERGPNTVHPSVAPISSEETPSPELVASGSQAPKIEPIHAT</p> <p>ADPQKLSVLITPVPDTQAATRVDHHHHHHH</p>
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only