

Recombinant Human Phosphofructokinase 1,PFKM,PFKX (C-6His)



Catalog No: #AP76018

Orders: order@signalwayantibody.com

Package Size: #AP76018-1 10ug #AP76018-2 50ug #AP76018-3 500ug #AP76018-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Phosphofructokinase 1,PFKM,PFKX (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P08237
Uniprot	P08237
GeneID	5213;
Calculated MW	86.1kDa
Target Sequence	THEEHHAAKTLGIGKAIIVLTSGGDAQGMNAAVRAVVRVVGIFTGARVFFVHEGYQGLVDGGDHIKEATWESV SMMLQLGGTVIGSARCKDFRERERGLRAAYNLVKGITNLCVIGGDGSLTGADTFRSEWSDLLSDLQKAGKIT DEEATKSSYLNIVGLVGSIDNDFCGTDMTIGTDSALHRIMEIVDAITTTAQSHQRTFVLEVMMGRHCGYLALVTSL SCGADWVFIPECPPDDDDWEEHLCRRLSETRTRGSRLNIIIVAEGAIKNGKPIITSEDIKNLVVKRLGYDTRVTVL GHVQRGGTPSAFDRILGSRMGVEAVMALLEGTPDTPACVVSLSGNQAVRLPLMECVQVTKDVTKAMDEKKF DEALKLRGRSFMNNWEVYKLLAHRPPVSKSGSHTVAVMNVGAPAAAGMNAAVRSTVIRIGLIQGNRVLVVDH GFEGGLAKGQIEEAGWSYVGGWGTGQGGSKLGTKRTLPKKSFEQISANITKFNIQGLVIIGGFEAYTGGLELMEG RKQFDELCPFFVIPATVSNNVPGSDFSVGADTALNTICTTCDRIKQSAAGTKRRRVIIETMGGYCGYLATMAGL AAGADAAYIFEEPFTIRDLANVEHLVQKMKTTVKRGLVLRNEKCNENYTTDFIFINLYSEEKGFIDSRKNVGL HMQQGGSPFPDRNFATKMGAKAMNWMMSGKIKESYRNGRIFANTPDSGCVLGMRKRALVFQPVAVELKDQT DFEHRIPKEQWWLKLRLPILKILAKYEIDLDTSDHAHLEHITRKRSGEAAVVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only