

Recombinant Human Fructose-2,6-Bisphosphatase 1,PFKFB1,PFK,FBPase 1 (C-6His)



Catalog No: #AP76020

Orders: order@signalwayantibody.com

Package Size: #AP76020-1 10ug #AP76020-2 50ug #AP76020-3 500ug #AP76020-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Fructose-2,6-Bisphosphatase 1,PFKFB1,PFK,FBPase 1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P16118
Uniprot	P16118
GeneID	5207;
Calculated MW	55.6kDa
Target Sequence	SPEMGELTQTRLQKIWIPIHSSGSSRLQRRRGSSIPQFTNSPTMVMVGLPARGKTYISTKLTRYLNWIGTPTKV FNLGQYRREAVSYKNEYEFFLPDNMEALQIRKQCALAALKDVHNYLSHEEGHVAVFDATNTTRRRSLILQFAK EHGYKVFIESICNDPGIIAENIRQVKLGSPDYIDCDREKVLDFLKRIECYEVNYQPLDEELDSHLSYIKIFDVGT RYMVNRVQDHIQSRTVYYLMNIHVTPRSIYLGRHGESELNIRGRIGGDSGLSVRGKQYAYALANFIQSQGISL KVWTSHMKRTIQTAELGVPYEQWKALNEIDAGVCEEMTYEEIQEHYPEEFALRDQDKYRYRYPKGESYEDL VQRLEPVIMELERQENVLVICHQAVMRCLLAYFLDKSSDELPLYKCPHHTVTKLTPVAYGCKVESIYLNVEAVNT HREKPEVDITREPEEALDTPAHYVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only