

Recombinant Human Phosphoglycerate Kinase 1,PGK1 (C-6His)

Catalog No: #AP76026

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76026-1 10ug #AP76026-2 50ug #AP76026-3 500ug #AP76026-4 1mg

Description

Product Name	Recombinant Human Phosphoglycerate Kinase 1,PGK1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P00558
Uniprot	P00558
GeneID	5230;
Calculated MW	45.5kDa
Target Sequence	SLSNKLTLDKLDVKGKRVVMRVDNFNPMKNNQITNNQRIKAAVPSIKFCLDNGAKSVVLMShLGRPDGVPMP DKYSLEPVAVELKSLGKDVFLKDCVGPEVEKACANPAAGSVILLENLRFHVVEEGKGDASGNKVKAEPAKI EAFRASLSKLGDVYVNDAFGTAHRAHSSMVGVNLPQKAGGFLMKKELNYFAKALESPPERPFLAILGGAKVAD KIQLINMLDKVNEMIIGGGMAFTFLKVLNNMEIGTSLFDEEGAKIVKDLMSKAEKNGVKITLPVDFVTADKFDE NAKTGQATVASGIPAGWMGLDCGPESKKYAEAVTRAKQIVWNGPVGVFWEAFARGTKALMDEVVKATSR GCITIIGGGDTATCCAKWNTEDKVSHVSTGGGASLELLEGKVLPGVDALSNIIDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris,150mM NaCl,20% glycerolpH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only