

Recombinant Human IL-17 Receptor A,IL-17RA,CD217 (C-Fc)

Catalog No: #AP76027

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76027-1 10ug #AP76027-2 50ug #AP76027-3 500ug #AP76027-4 1mg

Description

Product Name	Recombinant Human IL-17 Receptor A,IL-17RA,CD217 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q96F46
Uniprot	Q96F46
GeneID	23765;
Calculated MW	60.6kDa
Target Sequence	LRLLDHRALVCSQPGLNCTVKNSTCLDDSWIHPRNLTPSSPKDLQIQLHFAHTQQGDLFPVAHIEWTLQTDASILYLEGAELSVLQLNTNERLCVRFEFLSKLRHHHRRWRFTFSHFVVDPDQEYEVTVHHLPKPIPDGDPNHQSKNFLVPDCEHARMKVTTPCMSSGSLWDPNITVETLEAHQLRVSF TLWNETHYQILLTSFPHMENHSCFEHMHHPAPRPEEFHQRSNVTLTLRNLKGCCRHQVQIQPFSSCLNDCLRHSATVSCPEMPDTPPEPIPDYMLPWDDIEGRMDEPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNVFSCSV MHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only