

Recombinant Human Cystatin C,CST3 (C-6His)

Catalog No: #AP76029



Package Size: #AP76029-1 10ug #AP76029-2 50ug #AP76029-3 500ug #AP76029-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Cystatin C,CST3 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P01034
Uniprot	P01034
GeneID	1471;
Calculated MW	14.4kDa
Target Sequence	SSPGKPPRLVGGPMDASVEEEGVRRALDFAVGEYNKASNDMPYHSRALQVVRARKQIVAGVNYFLDVELGRT TCTKTQPNLDNCPFHDQPHLKRKAFCSFQIYAVPWQGTMTLSKSTCQDAVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 10mM PB, 200mM Nacl, pH6.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only