

## Recombinant Human Interleukin-5,IL-5 (C-6His)

Catalog No: #AP76031



Package Size: #AP76031-1 10ug #AP76031-2 50ug #AP76031-3 500ug #AP76031-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Interleukin-5,IL-5 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P05113
Uniprot	P05113
GeneID	3567;
Calculated MW	14.2kDa
Target Sequence	IPTEIPTSALVKETLALLSTHRTLIIANETLRIPVPVHKNHQLCTEEIFQGGIGTLESQTVQGGTVERLFKNLSLIKKY IDGQKKKCGEERRRVNQFLDYLQEFLGVMNTEWIIIESVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only