

Recombinant Human Thymidine Kinase, Cytosolic,TK1 (C-6His)

Catalog No: #AP76035

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76035-1 10ug #AP76035-2 50ug #AP76035-3 500ug #AP76035-4 1mg

Description

Product Name	Recombinant Human Thymidine Kinase, Cytosolic,TK1 (C-6His)
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P04183
Uniprot	P04183
GeneID	7083;
Calculated MW	26.5kDa
Target Sequence	MSCINLPTVLPGSPKTRGQIQVILGPMFSGKSTELMRRVRRFQIAQYKCLVIKYAKDTRYSSSFCTHDRNTME ALPACLLRDVAQEALGVAVIGIDEGQFFPDIVFCEAMANAGKTVIVAALDGTFRKPFGAILNLVPLAESVVKL TAVCMECFREAAAYTKRLGTEKEVEVIGGADKYHSVCRLCYFKKASGQPAGPDNKENCVPVPGKPGEAVAARKL FAPQQILQCSPANVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of PBS,pH7.4.
Storage	Store at -70C.

Note: This product is for in vitro research use only