

Recombinant Mouse CD28,,TP44 (C-Fc-6His)

Catalog No: #AP76037



Package Size: #AP76037-1 10ug #AP76037-2 50ug #AP76037-3 500ug #AP76037-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse CD28,,TP44 (C-Fc-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P31041
Uniprot	P31041
GeneID	12487;
Calculated MW	43kDa
Target Sequence	NKILVKQSPLLVDVDSNEVSLSCRYSYNLLAKEFRASLYKGVNSDVEVCVGNNGNFTYQPQFRSNAEFNCDGDF DNETVTFRLWNLHVNHTDIYFCKIEFMYPPPYLDNERSNGTIIHIKEKHLCHTQSSPKVDDIEGRMDEPKSCDK THTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREE QYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCL LVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNVFNCSVMHEALHNHYTQ KSLSLSPGKHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only