

Recombinant Mouse 4-1BB,TNFRSF9,CD137 (C-Fc)

Catalog No: #AP76039



Package Size: #AP76039-1 10ug #AP76039-2 50ug #AP76039-3 500ug #AP76039-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse 4-1BB,TNFRSF9,CD137 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P20334
Uniprot	P20334
GeneID	21942;
Calculated MW	44.7kDa
Target Sequence	VQNSCDNCQPGTFCRKYNPVCKSCPPSTFSSIGGQPNCNICRVCAGYFRFKKFCSSSTHNAECECIEGFHCLG PQCTRCEKDCRPGQELTKQGCKTCSLGTFNQNGTGVCRPWTNCSLDGRSVLKTGTTEKDVCGPPVVSF SPSTTISVTPEGGPGGHSLLQVLVDDIEGRMDEPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPE EVTQCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKAL PAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSD GSFFLYSKLTVDKSRWQQGNVFCSCVMHEALHNHYTQKSLSLSPGK
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only