

# Recombinant Human Sedoheptulokinase, SHPK (C-6His)

Catalog No: #AP76056

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP76056-1 10ug #AP76056-2 50ug #AP76056-3 500ug #AP76056-4 1mg

## Description

Product Name	Recombinant Human Sedoheptulokinase, SHPK (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9UHH6
Uniprot	Q9UHH6
GeneID	23729;
Calculated MW	52.5kDa
Target Sequence	MAARPITLGIDLGTTSVKAALLRAAPDDPSGFVAVLASCARAARAEAAVESAVAGPQGREQDVSRIQLQALHECLA ALPRPQLRSVVGIGVSGQMHGVMVFWKTGQGCWEWTEGGITPVFEPRAVSHLVTWQDGRCSSEFLASLPQPKS HLSVATGFGCATIFWLLKYRPEFLKSYDAAGTIHDYVAMLCGLPRPLMSDQNAASWGYFNTQSQSWNVETL RSSGFVHLLPDIAEPGSVAGRTSHMWFEIPKGTQVGVALGDLQASVYSCMAQRDVAVLNISTSVQLAASMP SGFQPAQTPDPTAPVAYFPYFNRTYLGAASLNGGNVLAFLVHMLVQWMADLGLVEESTVYSRMIQAAVQ QRDTHLTITPTVLGERHLPDQLASVTRISSDLSLGHVTRALCRGIVQNLHSMPLPIQLQEWGVERVMGSGSA LSRNDVLKQEVQRAFPLPMSFGQDVDAVGAALVMLRRHLNPKESVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only