

Recombinant Human Selenophosphate Synthase 1,SEPHS1 (C-6His)

Catalog No: #AP76057

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76057-1 10ug #AP76057-2 50ug #AP76057-3 500ug #AP76057-4 1mg

Description

Product Name	Recombinant Human Selenophosphate Synthase 1,SEPHS1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P49903
Uniprot	P49903
GeneID	22929;
Calculated MW	43.9kDa
Target Sequence	MSTRESFNPESEYELDKSFRLTRFTELKGTGCKVPQDVLQKLLLESLQENHFQEDEQFLGAVMPRLGIGMDTCVI PLRHGGLSLVQTTDIYPIVDDPYMMGRIACANVLSLDLYAMGVTECDNMLMLLGVSNKMTDRERDKVMPLIQ GFKDAAEEAGTSVTGGQTVLNPWIVLGGVATTVCQPNEFIMPDNAVPGDVLVLTkPLGTQVAVAVHQWLDIP EKWNKIKLVVTQEDVELAYQEAMMNMARLNRTAAGLMHTFNAAHAATDITGFGILGHAQNLAQQRNEVSFVIH NLPVLAKMAAVSKACGNMFMHGTCPETSGLLICLPREQAARFCAEIKSPKYGEGHQAWIIGIVEKGNRTA RIIDKPRIIEVAPQVATQNVNPTPGATSVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 25mM Tris.HCl,pH7.3,100mM glycine,10% glycerol.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only