

Recombinant Human Receptor Tyrosine Kinase MerTK,MERTK,MER (C-6His)



Catalog No: #AP76058

Orders: order@signalwayantibody.com

Package Size: #AP76058-1 10ug #AP76058-2 50ug #AP76058-3 500ug #AP76058-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Receptor Tyrosine Kinase MerTK,MERTK,MER (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q12866
Uniprot	Q12866
GeneID	10461;
Calculated MW	36kDa
Target Sequence	MKINNEEIVSDPIYIEVQGLPHFTKQPESMNVTRNTAFNLTCQAVGPPEPVNIFVWQNSSRVNEQPEKSPSVLT VPGLTEMAVFSCEAHNDKGLTVSKGVQINIKAIKIPSPTEVSIRNSTAHSILISWVPGFDGYSPFRNCSIQVKEAD PLSNGSVMIFNTSALPHLYQIKQLQALANYSIGVSCMNEIGWSAVSPWILASTTEGAPSVAPLNVTVFLNESSD NVDIRWMKPPTKQQDGELVGYRISHVWQSAGISKELLEEVGQNGSRARISVQVHNATCTVRIAAVTKGGVGP FSDPVKIFIPAHGWVDYAPSSTPAPGNAVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, pH8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only