

# Recombinant Human Coagulation Factor IX,F9 (C-6His)

Catalog No: #AP76059

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP76059-1 10ug #AP76059-2 50ug #AP76059-3 500ug #AP76059-4 1mg

## Description

Product Name	Recombinant Human Coagulation Factor IX,F9 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P00740
Uniprot	P00740
GeneID	2158;
Calculated MW	49.8kDa
Target Sequence	TVFLDHENANKILNRPKRYNSGKLEEFVQGNLERECMEEKCSFEEAREVFENTERTEFWKQYVDGDQCESN PCLNGGSKDDINSYECWCPFGFEGKNCELDVTCKNIKNGRCEQFCCKNSADNKVVCSCTEGYRLAENQKSCE PAVPFPCGRVSVSQTSLTRAETVFPDQVDYVNSTEAETILDNITQSTQSFNDFTRVVGEDAKPGQFPWQVVL NGKVDAFCGGSIVNEKWIVTAAHCVETGVKITVAGEHNIEETEHEQKRNIRIIPHHNYNAAINKYNHDIALLE LDEPLVLSYVTPICIAADKEYTNIFLKFSGYVSGWGRVVFHKGRSALVLQYLRVPLVDRATCLRSTKFTIYNMFM CAGFHEGGRDSCQGDSSGPHVTEVEGTSFLTGIISWGEECAMKGKYGITKVSRYVNWIKETKLTVDHHHH HH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl,150mM NaCl,10% Glycerol,pH8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only