

Recombinant Mouse BMP Receptor IA,ALK-3,CD292 (C-Fc-6His)



Catalog No: #AP76066

Orders: order@signalwayantibody.com

Package Size: #AP76066-1 10ug #AP76066-2 50ug #AP76066-3 500ug #AP76066-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse BMP Receptor IA,ALK-3,CD292 (C-Fc-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P36895
Uniprot	P36895
GeneID	12166;
Calculated MW	42.2kDa
Target Sequence	QNLDSMLHGTGMKSDLQKKPENGVTLAPEDTLFPLKCYCSGHCPDDAINNTCITNGHCFIIEEDDQGETTL TSGCMKYEGSDFQCKDSPKAQLRRRTIECCRTNLCNQYLQPTLPPVIGPFFDGSIRVDDIEGRMDEPKSCDKT HTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQ YNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCL VKGFYPSDIAVEWESNGQPENNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNVFCFSVMHEALHNHYTQK SLSLSPGKHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only