

# Recombinant Mouse Plasma Glutamate Carboxypeptidase, PGCP (C-6His)



Catalog No: #AP76072

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP76072-1 10ug #AP76072-2 50ug #AP76072-3 500ug #AP76072-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Mouse Plasma Glutamate Carboxypeptidase, PGCP (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q9WVJ3
Uniprot	Q9WVJ3
GeneID	54381;
Calculated MW	50.9kDa
Target Sequence	KAVFKNGVSQRTFREIKEEIANYEDVAKAIINLAVYGKYQNRSYERLGLLVDTVGPRLSGSKNLEKAIQIMYQNL QQDGLNVHLEQVRIPHWERGEESAVMLEPRIHKMAILGLGSSIGTPPGGITAEVLVASFDELQRRASEARG KIIVYNQPYTYEKTQYRVQGAVEAAKVGAVASLIQSVASFISIYSPHTGIQKYQDGVPKIPTACITVEDAEMMS RMSRGNKIVIHLEMGAKTYPDTDSFNTVAEITGSMYPPEEVLVSGHLDSDWDVGGALDDGGGAFISWEALS LVKDLGLRPKRTLRLVLWTAEEQGGIGASQYYELHKANISKYSLVMEADSGTFLPTGLQFTGSDKARAIMKEV MNLQPLNVTKVFSNGEGTDINFWIQAGVPGASLRDDLKYFFFHSHGDTMTVMDPKQMNVAADVAVVA YVVADMDEMLPRSDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only