

Recombinant Human Cadherin-16,CDH16 (C-6His)

Catalog No: #AP76077



Package Size: #AP76077-1 10ug #AP76077-2 50ug #AP76077-3 500ug #AP76077-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Cadherin-16,CDH16 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O75309
Uniprot	O75309
GeneID	1014;
Calculated MW	84.6kDa
Target Sequence	MPKAQPAELSVEVPENYGGNFPLYLTKLPLPREGAEGQIVLSGDSGKATEGPFAMDPDPSGFLLVTRALDREE QAEYQLQVTLEMQDGHVWLPVQVPLVHVKDENDQVPHFSQAIYRARLSRGRTRPGIPFLFLEASDRDEPGTAN SDLRFHILSQAPAPSPDMFQLEPRLGALALSPKGSTSLDHALERTYQLLVQVKDMGDQASGHQATATVEVSI IESTWWSLEPIHLAENLKVLYPHHMAQVHWSGGDVHYHLESHPPGPFVNAEENLYVTRELDREAQAEYLLQ VRAQNSHGEDIYAAPLELHVLVMDENDNVPICPPRDPTVSIPELSPPGTEVTRLSAEDADAPGSPNSHVYQLL SPEPEDGVEGRAQVDPTSGSVTLGVLPLRAGQNILLVLAMDLAGAEGGFSSTCEVEVAVTDINDHAPEFITS QIGPISLPEDVEPGTLVAMLTADADLEPAFRLMDFAIERGDTEGTFGLDWEPPDSGHVRLRLCKNLSYEAAPSH EVVVVVQSVAKLVGPGPGGATATVTVLVERVMPPPDKDQESYEASVPISAPAGSFLTIQPSDPISRTLRFSL VNDSEGWLCIEKFSGEVHTAQSLQGAQPGDITYVLEAQDTDEPRLSASAPLVIHFLKAPPAPALTAPVPSQ YLCTPRQDHGLIVSGPSKDPDLASGHGYPYSFTLGNPTVQRDWRLQTLNGSHAYLTLALHWVEPREHIIPVVV SHNAQMWQLLVRVIVCRCNVEGQCMRKVGRMKGMPTKLSAVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only