

# Recombinant Human Prion-Like Protein Doppel,PRND (C-6His)

Catalog No: #AP76082

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP76082-1 10ug #AP76082-2 50ug #AP76082-3 500ug #AP76082-4 1mg

## Description

Product Name	Recombinant Human Prion-Like Protein Doppel,PRND (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9UKY0
Uniprot	Q9UKY0
GeneID	23627;
Calculated MW	15.5kDa
Target Sequence	RGIKHRIKWNRKALPSTAQITEAQVAENRPGAFIKQGRKLDIDFGAEGNRYYEANYWQFPDGIHYNGCSEANV TKEAFVTGCINATQAANQGEFQKPDNKLHQVWLWRLVQELCSLKHCEFWLERGVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only