

# Recombinant Human Interleukin-12 Subunit $\beta$ ,IL-12 p40,IL-12B



Catalog No: #AP76085

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP76085-1 10ug #AP76085-2 50ug #AP76085-3 500ug #AP76085-4 1mg

## Description

Product Name	Recombinant Human Interleukin-12 Subunit $\beta$ ,IL-12 p40,IL-12B
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P29460
Uniprot	P29460
GeneID	3593;
Calculated MW	34.6kDa
Target Sequence	IWELKKDVYVVELDWYPDAPGEMVVLTCDTPEEDGITWTLTLDQSSEVLGSGKTLTIQVKEFGDAGQYTCHKGG EVLSHSLLLLHKKEDGIWSTDILKDQKEPKNKTFLRCEAKNYSGRFTCWLLTISTDLTFSVKSSRGSSDPQG VTCGAATLSAERVRGDNKEYEYSVEQCEDSACPAAEESLPIEVMVDAVHKLKYENYTSSFFIRDIKPDPPKNL QLKPLKNSRQVEVSWEYPDTWSTPHSYFSLTFCVQVQGKSKREKKDRVFTDKTSATVICRKNASISVRAQDR YYSSWSEWASVPCS
Formulation	Lyophilized from a 0.2 $\mu$ m filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only