

Recombinant Human Kelch-Like protein 41,Kell,CD238 (C-6His)

Catalog No: #AP76086

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76086-1 10ug #AP76086-2 50ug #AP76086-3 500ug #AP76086-4 1mg

Description

Product Name	Recombinant Human Kelch-Like protein 41,Kell,CD238 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P23276
Uniprot	P23276
GeneID	3792;
Calculated MW	78.4kDa
Target Sequence	NFQNCGPRPCETSVCLDLRDHYLASGNTSVAPCTDFFSFACGRAKETNNSFQELATKNKNLRRRILEVQNSW HPGSGEEKAFQFYNSCMDTLAIEAAGTGPLRQVIEELGGWRISGKWTSLNFNRTLRLLSQYGHFPFFFRAYL GPHPASPHTPVIQIDQPEFDVPLKQDQEQKIYAQIFREYLTYLNQLGTLGGDPSKVQEHSSLSISITSRLFQFL RPLEQRRRAQGKLFQMVTIDQLKEMAPIDWLSCLQATFTPMSPSQSLVVDVEYLKNMSQLVEEMLLKQR DFLQSHMILGLVVTLSPALDSQFQEARRKLSQLRELTEQPPMPARPRWMKCVEETGTFEPTLAALFVREAF GPSTRSAAMKLFTAIRDALITRLRNLPMNEETQNMAQDKVAQLQVEMGASEWALKPELARQEYNDIQLGSS FLQSVLSCVRSRLRARIVQSFLQPHQHRWKVSPWDVNAYYSVSDHVVVFPAGLLQPPFFHPGYPRAVNFCAA GSIMAHHELLHIFYQLLLPGGCLACDNHALQEAHLCLKRHYAAFPLPSRTSFNDSLTFLENAADVGGGLAIALQAY SKRLLRHHGETVLPDLSPQIIFRSYAQVMCRKPSQDSDHDSHPPHLRVHGPLSSTPAFARYFRCARGA LLNPSSRCQLWVDHHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only