

Recombinant Human Serpin B9, SERPINB9 (C-6His)

Catalog No: #AP76088



Package Size: #AP76088-1 10ug #AP76088-2 50ug #AP76088-3 500ug #AP76088-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Serpin B9, SERPINB9 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P50453
Uniprot	P50453
GeneID	5272;
Calculated MW	43.4kDa
Target Sequence	METLSNASGTFAIRLLKILCQDNPSHNVFCSPVSISSALAMVLLGAKGNTATQMAQALSLNTEEDIHRAFQSLT EVNKAGTQYLLRTANRLFGEKTCQFLSTFKESCLQFYHAELKELSFIRAAEESRKHINTWVSKKTEGKIEELLP GSSIDAETRLVLVNAIYFKGKWNEPFDETYREMPFKINQEEQRPVQMMYQEATFKLAHVGEVRAQLLELPYA RKELSLVLLPDDGVELSTVEKSLTFEKLTAWTKPDCMKSTEVEVLLPKFKLQEDYDMESVLRHLGIVDAFQQ GKADLSAMSAERDLCLSKFVHKSFVEVNEEGTEAAAASSCFVVAECCMESGPRFCADHPFLFFIRHNRANSIL FCGRFSSPVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only