

Recombinant Human Interleukin-23,IL-23 (C-6His)

Catalog No: #AP76092



Package Size: #AP76092-1 10ug #AP76092-2 50ug #AP76092-3 500ug #AP76092-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Interleukin-23,IL-23 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9NPF7 & P29460
Uniprot	Q9NPF7
GeneID	51561;
Calculated MW	54.4kDa
Target Sequence	RAVPGGSSPAWTQCQQLSQKLCTLAWSAHLVGHMDLREEGDEETTNDVPHIQCGDGCDDPQGLRDNSQFC LQRIHQGLIFYEKLLGSDIFTGEPSELLPDSVPGQLHASLLGLSPLLQPEGHHWETQQIPSLSPSQPWQRLLLR KILRSLQAFVAVAAARVFAHGAATLSPVDHHHHH&IWELKDKDVYVVELDWYPDAPGEMVVLTCDTPEEDGITW TLDQSSEVLGSGKTLTIQVKEFGDAGQYTCHKGGEVLSHLLLLHKKEDGIWSTDILKDQKEPKNKTFRLRCEAK NYSGRFTCWLLTTISTDLTFSVKSSRGSSDPQGVTCGAATLSAERVRGDNKEYEYSVEQCEDSACPAAEESL PIEVMVDAVHKLKYENYSSFFIRDIIKPDPPKNLQLKPLKNSRQVEVSWEYPTDWSTPHSYFSLTFCVQVQGGK SKREKKDRVFTDKTSATVICRKNASISVRAQDRYSSSWSEWASVPCS
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only