

Recombinant Mouse Tyrosine-Protein Kinase Receptor TYRO3,Dtk (C-mFc)



Catalog No: #AP76096

Orders: order@signalwayantibody.com

Package Size: #AP76096-1 10ug #AP76096-2 50ug #AP76096-3 500ug #AP76096-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Tyrosine-Protein Kinase Receptor TYRO3,Dtk (C-mFc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P55144
Uniprot	P55144
GeneID	22174;
Calculated MW	68.2kDa
Target Sequence	AGLKLMGAPVKMTVSQQQPVKLNCSVEGMEDPDIHWMKDGTVVQNASQVSISEHSWIGLLSLKSVERSDA GLYWCQVKDGEETKISQSVWLTVEGVPFFTVEPKDLAVPPNAPFQLSCEAVGPPPEPTIYWWRGLTKVGGPA PSPSVLNVTGVTQRTEFSCEARNIKGLATSRPAIVRLQAPPAAPFNNTVTTISSYNASVAWVPGADGLALLHSC TVQVAHAPGEWEALAVVVPVPPFTCLLRNLAPATNYSLRVRCANALGSPSYGDWVPFQTKGLAPARAPQNF HAIRTDGLILEWEEVIPEDPGEGPLGPYKLSWVQENGTQDELMVEGTRANLTDWDPQKDLILRVCASNAIGD GPWSQPLVSSHDHAGRQGGPPHSRTSGGGGSGGGGSGGGGSPRDCGCKPCICTVPEVSSVFIFPPKPKDV LTITLTPKVTCVVDISKDDPEVQFSWFVDDVEVHTAQTQPREEQFNSTFRSVSELPIMHQDWLNGKEFKCRV NSAAFPAPIEKTISKTKGRPKAPQVYTIPPPKEQMAKDKVSLTCMITDFFPEDITVEWQWNGQPAENYKNTQPI MDTDGSYFVYSKLNQKSNWEAGNTFTCSVLHEGLHNHHTEKSLSHSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only