

Recombinant Human Resistin,RETN,ADSF (C-Fc)

Catalog No: #AP76097



Package Size: #AP76097-1 10ug #AP76097-2 50ug #AP76097-3 500ug #AP76097-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Resistin,RETN,ADSF (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9HD89
Uniprot	Q9HD89
GeneID	56729;
Calculated MW	36.8kDa
Target Sequence	<p>KTLCSMEEAINERIQEVAGSLIFRAISSIGLECSVTSRGLATCPRGFAVTGCTCGSACGSWDVRAETTCHC</p> <p>QCAGMDWTGARCCRVQPVDIEGRMDEPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTC</p> <p>VVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIE</p> <p>KTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTPPVLDSDGSFFL</p> <p>YSKLTVDKSRWQQGNVFCSSVMHEALHNHYTQKSLSLSPGK</p>
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only