

Recombinant Human IL-2 Receptor Subunit γ ,IL-2RG,CD132 (C-Fc)

Catalog No: #AP76098

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76098-1 10ug #AP76098-2 50ug #AP76098-3 500ug #AP76098-4 1mg

Description

Product Name	Recombinant Human IL-2 Receptor Subunit γ ,IL-2RG,CD132 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P31785
Uniprot	P31785
GeneID	3561;
Calculated MW	55.4kDa
Target Sequence	LNTILTTPNGNEDTTADFFLTTMPTDSLVSSTLPLPEVQCFVFNVEYMNCTWNSSSEPQPTNLTLYWYKNSD NDKVQKCSHYLFSEEITSGCQLQKKEIHLYQTFVVLQDPPREPRRQATQMLKLQNLVIPWAPENLTLHLKLS QLELNWNNRFLNHCLEHLVQYRTDWDHSWTEQSVDYRHKFSLPSVDGQKRYTFRVRSRFPNPLCGSAQHWS EWSHPHGWGNTSKENPFLFALEAVDDIEGRMDEPKSCDKTHTCPPCPAPELGGPSVFLFPPKPKDTLMISR TPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNK ALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVLD SDGSFFLYSKLTVDKSRWQQGNVFCFSVMHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 μ m filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only