

Recombinant Human Lymphocyte Antigen 86,LY86,MD-1 (C-Fc)

Catalog No: #AP76100

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76100-1 10ug #AP76100-2 50ug #AP76100-3 500ug #AP76100-4 1mg

Description

Product Name	Recombinant Human Lymphocyte Antigen 86,LY86,MD-1 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O95711
Uniprot	O95711
GeneID	9450;
Calculated MW	42.8kDa
Target Sequence	GGGGKAWPTHVVCSDSGLEVLYQSCDPLQDFGFSVEKCSKQLKSNIRFGIILREDIKELFLDLALMSQGSSV LNFSYPICEAALPKFSFCGRRKGEQIYYAGPVNNPEFTIPQGEYQVLELYTEKRSTVACANATIMCSVDIEG RMDEPKSCDKTHTCPPCPAPPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVDVSHEDPEVKFNWYVDGVEV HNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREE MTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSGDSFFLYSKLTVDKSRWQQGNVFCSSVM HEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only