

Recombinant Human Aldehyde Dehydrogenase 3-A1,ALDH3A1 (C-6His)

Catalog No: #AP76103

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76103-1 10ug #AP76103-2 50ug #AP76103-3 500ug #AP76103-4 1mg

Description

Product Name	Recombinant Human Aldehyde Dehydrogenase 3-A1,ALDH3A1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P30838
Uniprot	P30838
GeneID	218;
Calculated MW	51.4kDa
Target Sequence	MSKISEAVKRARAAFSSGRTRPLQFRIQQLEALQRLIQEQEQELVGALAADLHKNEWNAYEEVVYVLEEIEY MIQKLPEWAADPEVEKTPQTQQDELYIHSEPLGVVLVIGTWNYPFNLTIQPMVGAIAAGNAVVLKPSSESENMA SLLATIIPQYLDKDLYPVINGGVPETTELLKERFDHILYTGSTGVGKIIMTAAAKHLTPVTLELGGKSPCYVDKNC DLDVACRRIAWGKFMNSGQTCVAPDYILCDPSIQNQIVEKLLKSLKEFYGEDAKKSRDYGRIISARHFQRMGL IEGQKVAYGGTGDAATRYIAPTILTDVDPQSPVMQEEIFGPVLPVIVCVRSLLEEAIQFINQREKPLALYMFSSNDK VIKKMIAETSSGGVAANDVIVHITLHSLPFGGNGSGMGSYHGKKSFEFESHRRSCLVRPLMNDEGLKVRYP SPAKMTQHVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only