

# Recombinant Human IL-1 Receptor-Like 2,IL-1RL2 (C-Fc)



Catalog No: #AP76104

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP76104-1 10ug #AP76104-2 50ug #AP76104-3 500ug #AP76104-4 1mg

## Description

Product Name	Recombinant Human IL-1 Receptor-Like 2,IL-1RL2 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9HB29
Uniprot	Q9HB29
GeneID	8808;
Calculated MW	63.1kDa
Target Sequence	DGCKDIFMKNEILSASQPFAFNCTFPPITSGEVSVTWYKNSSKIPVSKIIQSRIHQDETWILFLPMEWGDSDGVYQ CVIKGRDSCRIHVNLTVFEKHWCDTSIGGLPNLSDEYKQILHLGKDDSLTCHLHFPKSCVLGPIKWKDCNEI KGERFTVLETRLLVSNVSAEDRGNYACQAILTHSGKQYEVLNGITVSITERAGYGGSVPKIIPKNHSIEVQLGT TLIVDCNVTDTKDNTNLRCWRVNNTLVDDYDESKRIREGVETHVSVFREHNLYTVNITFLEVKMEDYGLPFMC HAGVSTAYIILQLPAPDFRAYDIEGRMDEPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTC VVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIE KTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTPPVLDSDGSFFL YSKLTVDKSRWQQGNVFCFSVMHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only