

Recombinant Mouse Semaphorin 4C,SEMA4C (C-Fc)

Catalog No: #AP76108



Package Size: #AP76108-1 10ug #AP76108-2 50ug #AP76108-3 500ug #AP76108-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Semaphorin 4C,SEMA4C (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q64151
Uniprot	Q64151
GeneID	20353;
Calculated MW	99.3kDa
Target Sequence	AEMWWNLVPRKTVSSGELVTVRRRFSQTGIQDFLTLTEHSGLLYVGAREALFAFSVEALELQGAISWEAPA EKKIECTQKGSNQTECFNIRFLQPYNSSHLVCGTYAFQPKCTYINMLTFTLDRAEFEDGKKGKCPYDPAKG HTGLLDVDELYSATLNNFLGTPEVILRYMGTHHSIKTEYLAFWLNEPHFVGSFVPEVSGSFTGDDDKIYFFFS ERAVEYDCYSEQVVARVARVCKGDMGGARTLQKKWTTFLKARLVCSAPDWKVFYFNQLKAVHLRGASWHN TTFGFVFQARWGMDLSAVCEYQLEIQQVFEGPYKEYSEQAQKWARYTDPVSPRPGSCINNWHRDNGY TSSLELPDNTLNFIKKHPLMEDQVKPRLGRPLLVKKNNTNFTHVADRVPGLDGATYTVLFIGTGDGWLLKAVSL GPWIHMVEELQVFDQEPVESLVSQSCKVLFAGSRSQLVQLSLADCTKYRFCVDCVLARDPYCAWNVNTRC VATTSGRSGSFLVQHVANLDTSKMCNQYGIKKVRISIPKNITVSGTDLVLPCHLSSNLAHAHWTFGSQDLPAE QPGSFLYDTGLQALVMAAQRHSGPYRCYSEEQGTRLAESAESYLVAVVAGSSVTLEARAPLENLGVDIEGR MDEPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVH NAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEM TKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTPPVLDSDGSFFLYSKLTVDKSRWQQGNVFCFSVMH EALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only