

# Recombinant Human Ephrin A Receptor 4,EphA4 (C-Fc)

Catalog No: #AP76112

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP76112-1 10ug #AP76112-2 50ug #AP76112-3 500ug #AP76112-4 1mg

## Description

Product Name	Recombinant Human Ephrin A Receptor 4,EphA4 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P54764
Uniprot	P54764
GeneID	2043;
Calculated MW	85.6kDa
Target Sequence	VTGSRVYPANEVTLLD SRSVQ GELGWIASPLEGGWEEVSIMDEKNTPIRTYQVCNVMEPSQNNWLRTDWITR EGAQRVYIEIKFTLRDCNSLPGVMGTCKETFNLYYESDNDKERFIRENQFVKIDTIAADESFTQVDIGDRIMKL NTEIRDVGPLSKKGFYLAQDVGACIALVSVRVFYKKCPLTVRNLAQFPDTITGADTSSLVEVRGSCVNNSEEK DVPKMYCGADGEWLVPIGNCLCNAGHEERSGECQACKIGYYKALSTDATCAKCPPHSYSVWEGATSCTCDR GFFRADNDAASMPCTRPPSAPLNLSNVNETSVNLEWSSPQNTGGRQDISYNNVCKKCGAGDPSKCRPCGS GVHYTPQQNGLKTTKVSITDLLAHTNYTFEIWAVNGVSKYNPNPDQSVSVTVTTNQAAPSSIALVQAKEVTRY SVALAWLEPDRPNGVILEYEVKYYEKDQNSRYRVRTAARNTDIKGLNPLTSYVFHVRARTAAGYGDFSEPL EVTNTVPSRIIGDGANSTVDDIEGRMDEPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTC VVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIE KTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSDGSFFL YSKLTVDKSRWQQGNVFCFSVMHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris,150mM NaCl,pH8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only