

Recombinant Human IL-2 Receptor Subunit α,IL-2RA,CD25 (C-Fc)

Catalog No: #AP76115

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76115-1 10ug #AP76115-2 50ug #AP76115-3 500ug #AP76115-4 1mg

Description

Product Name	Recombinant Human IL-2 Receptor Subunit α,IL-2RA,CD25 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P01589
Uniprot	P01589
GeneID	3559;
Calculated MW	48.9kDa
Target Sequence	ELCDDDPPEIPHATFKAMAYKEGTM LNCECKRGFRRIKSGSLYMLCTGNSSHSSWDNQCQCTSSATRNTTK QVTPQPEEQKERKTTEMQSPMQPVDQASLPGHCREPPWENEATERIYHFVVGQMVVYQCVQGYRALHRG PAESVCKMTHGKTRWTQPQLICTGEMETSQFPGEEKPQASPEGRPESETSCVDDIEGRMDEPKSCDKTHTC PPCAPPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNS TYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKG FYPSDIAVEWESNGQPENNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNV FSCSVMHEALHNHYTQKSLSL SPGK
Formulation	Lyophilized from a 0.2 μm filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only