

Recombinant Human Ephrin A Receptor 8,EphA8 (C-Fc)



Catalog No: #AP76116

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76116-1 10ug #AP76116-2 50ug #AP76116-3 500ug #AP76116-4 1mg

Description

Product Name	Recombinant Human Ephrin A Receptor 8,EphA8 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P29322
Uniprot	P29322
GeneID	2046;
Calculated MW	78.4kDa
Target Sequence	<p>ARGEVLLDSTIHGDWGWLTYPAHGWDSINEVDESFPQIHTYQVCNVMSPNQNNWLRTSWVPRDGARRV YAEIKFTLRDCNSMPGVLGTCKETFNLYLESDRDLGASTQESQFLKIDTIAADESFTGADLGVRRLKLNTEVR SVGPLSKRGFYLAQDQIGACLAILSLRIYYKKCPAMVRNLAASFSEAVTGADSSSLVEVRGQCVRHSEERDTPK MYCSAEGEWLVPIGKCVCSAGYEERRDACVACELGFYKSAPGDQLCARCPPHSHSAAPAAQACHCDLSYYR AALDPPSSACTRPPSAPVNLISVNGTSTLEWAPPLDPGGRSDITYNAVCRRCPPWALSRCACGSGTRFVP QQTSLVQASLLVANLLAHMNYSFWIEAVNGVSDLSPEPRRAAVNITTNQAGRRRNSVPQRPGPPASPASDP SRDQSSAGDVLWAFRQVPLWPCAPHQDPELEALHCLVDDIEGRMDEPKSCDKTHTCPPCPAPPELLGGPSVF LFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDW LNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQ PENNYKTTTPVLDSGSEFFLYSKLTVDKSRWQQGNVFCSSVMHEALHNHYTQKSLSLSPGK</p>
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only