

# Recombinant Mouse Legumain,Asparaginyl Endopeptidase (C-6His)

Catalog No: #AP76118

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP76118-1 10ug #AP76118-2 50ug #AP76118-3 500ug #AP76118-4 1mg

## Description

Product Name	Recombinant Mouse Legumain,Asparaginyl Endopeptidase (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	O89017
Uniprot	O89017
GeneID	19141;
Calculated MW	48.7kDa
Target Sequence	VPVGVDDPEDGGKHWVIVAGSNGWYNYRHQADACHAYQIIHRNGIPDEQIIVMMYDDIANSEENPTPGVVIN RPNGTDVYKGVLDYTGEDVTPENFLAVLRGDAEAVKKGSGKVLKSGPRDHVFIYFTDHGATGILVFPNDDL HVKDLNKTIRYMYEHKMYQKMFYIEACESGSMNHLRDDINVYATTAANPKESSYACYDEERGTYLGDWY SVNWMEDSDVEDLTKETLHKQYHLVKSHTNTSHVMQYGNKSISTMKVMQFQGMKHRASSPISLPPVTHLDLT PSPDVPLTILKRKLLRTNDVKESQNLIGQIQFLDARHVIEKSVHKIVSLLAGFGETAERHLSERTMLTAHDCYQ EAVTHFRTHCFNWHSVTYEHALRYLYVLANLCEAPYPIDRIEMAMDKVCLSHYVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris,150mM NaCl,20% Glycerol,pH8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only