

Recombinant Mouse Carbonic Anhydrase 14,CA14 (C-6His)

Catalog No: #AP76119

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76119-1 10ug #AP76119-2 50ug #AP76119-3 500ug #AP76119-4 1mg

Description

Product Name	Recombinant Mouse Carbonic Anhydrase 14,CA14 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q9WVT6
Uniprot	Q9WVT6
GeneID	23831;
Calculated MW	31.8kDa
Target Sequence	ADGGHHWTYEGPHGQDHWPTSYPECGGDAQSPINIQTDSVIFDPLPAVQPHGYDQLGTEPLDLHNNGHTV QLSLPPTLHLGGLPRKYTAAQLHLHWGQRGSLEGSEHQINSEATAELHVVHYDSQSYSSLSEAAQKPQGLA VLGILIEVGETENPAYDHILSRLEIRYKDQKTSVPPFSVRELFPPQLEQFFRYNGSLTTPPCYQSVLWTVFNR RAQISMGQLEKLQETLSSTEEDPSEPLVQNYRVPQLNQRTIFASFIQAGPLYTTGEMVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris,150mM NaCl,pH8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only