

## Recombinant Mouse PD-L1,B7-H1,CD274 (C-Fc)

Catalog No: #AP76123



Package Size: #AP76123-1 10ug #AP76123-2 50ug #AP76123-3 500ug #AP76123-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Mouse PD-L1,B7-H1,CD274 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q9EP73
Uniprot	Q9EP73
GeneID	60533;
Calculated MW	51.7kDa
Target Sequence	FTITAPKDLVVEYGSNVTMECRFPVERELDLLALVYWEKEDEQVIQFVAGEEDLKPQHSNFRGRASLPKDQ LLKGNAALQITDVKLQDAGVYCCIISYGGADYKRITLKVNPYRKINQRISVDPATSEHELICQAEGYPEAEVIWT NSDHQPVSGKRSVTTSRTEGMLLNVTSSLRVNATANDVFYCTFWRSQPGQNHTAELIPELPPATHPPQNRV DDIEGRMDEPKSCDKTHTCPPCPAPPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYV DGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLP PSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNVFS CSVMHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only