

Recombinant Human Lymphocyte Activation Gene 3,LAG-3,CD223 (C-6His)



Catalog No: #AP76125

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76125-1 10ug #AP76125-2 50ug #AP76125-3 500ug #AP76125-4 1mg

Description

Product Name	Recombinant Human Lymphocyte Activation Gene 3,LAG-3,CD223 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P18627
Uniprot	P18627
GeneID	3902;
Calculated MW	45.5kDa
Target Sequence	LQPGAEVPVWVAQEGAPAQLPCSPTIPLQDLSLLRRAGVTWQHQPDSGPPAAAPGHPLAPGPHPAAPSSWG PRPRRYTVLSVGGGLRSGRLPLQPRVQLDERGRQRGDFSLWLRPARRADAGEYRAAVHLRDRALSCRLRL RLGQASMTASPPGSLRASDWVILNCSFSRPDRPASVHWFRNRGQGRVPVRESPPHHHLAESFLFLPQVSPMD SGPWGCILTYRDGFNVSIMYNLTVLGLEPPTPLTVYAGAGSRVGLPCRLPAGVGTTRSFLTAKWTPPGGGPDL LVTGDNGDFTLRLEDVSAQAGTYTCHIHLQEQLNATVTLAIIITVTPKSFSGPSGLKLLCEVTPVSGQERFV WSSLDTPSQRSFSGPWLEAQEAQLLSQPWQCQLYQGERLLGAAVYFTELSSPGHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only