

Recombinant Human VEGF Receptor 2, VEGF R2, FLK-1, KDR (C-Fc)



Catalog No: #AP76126

Orders: order@signalwayantibody.com

Package Size: #AP76126-1 10ug #AP76126-2 50ug #AP76126-3 500ug #AP76126-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human VEGF Receptor 2, VEGF R2, FLK-1, KDR (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P35968
Uniprot	P35968
GeneID	3791;
Calculated MW	110.4kDa
Target Sequence	ASVGLPSVSLDLPRLSIQKDILTIKANTTLQITCRGQRDLDWLPNNQSGSEQRVEVTECSDGLFCKLTIPKVI GNDTGAYKCFYRETDLASVIYVYVQDYRSPFIASVSDQHGVVYITENKNKTVVIPCLGSISNLNVSLCARYPEK RFVPDGNRISWDSKKGFTIPSYMISYAGMVCFEAKINDESYQSIMYIVVVVGYRIYDVVLSPSHGIELSVGEKLV LNCTARTELVGIDFNWEYPSSKHQHKLVNRDLKTQSGSEMCKFLSTLTIDGITRSDQGLYTCAASSGLMTK KNSTFVRVHEKPFVAFGSGMESLVEATVGERVRIPAKYLGYPPEIKWYKNGIPLESNHTIKAGHVLTIMEVSE RDTGNYTVILTNPISKEKQSHVVSLVVYVPPQIGEKSLISPVDSYQYGTQTCTVYAIPPPHHIHWYWQLEEE CANEPHAVSVTNPYPCEEWRSVEDFQGGNKIEVNKNQFALIEGKNKTVSTLVIQAANVSALYKCEAVNKVGR GERVISFHVTRGPEITLQPDMPTEQESVSLWCTADRSTFENLTWYKLGQPPLPIHVGEPLTPVCKNLDLWK LNATMFSNSTNDILIMELKNASLQDQGDYVCLAQDRKTKKRHCVVRQLTVLERVAPTITGNLENQTTSIGESIEV SCTASGNPPPQIMWFKDNETLVEDSGIVLKDGNRNLTIKRVKEDGLYTCQACSVLGC AKVEAFFIIEGAQEK TNLEVDIEGRMDEPKCDKTHTCPPELGGPSVFLFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKF NWWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQ VYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTPPVLDSDGSFFLYSKLTVDKSRWQQ GNVFSCSVMHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only