

# Recombinant Human VEGF Receptor 1, VEGF R1, FLT-1 (C-Fc)

Catalog No: #AP76127

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP76127-1 10ug #AP76127-2 50ug #AP76127-3 500ug #AP76127-4 1mg

## Description

Product Name	Recombinant Human VEGF Receptor 1, VEGF R1, FLT-1 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P17948
Uniprot	P17948
GeneID	2321;
Calculated MW	109kDa
Target Sequence	SKLKDPELSLKGQTQHIMQAGQTLHLQCRGEEAAHKWSLPEMVSKESERLSITKSACGRNGKQFCSTLTLNNTAQ ANHTGFYSCKYLAVPTSKKKETESAIYIFISDTGRPFVEMYSEIPEIIHMTEGRELVIPCRVTSNPITVTLKKFPLD TLIPDGKRIIWDNRKGFIIISNATYKEIGLLTCEATVNGHLYKTNLTHRQTNTIIDVQISTPRPVKLLRGHTLVLNC TATPLNTRVQMTWSYPDEKNKRASVRRRIDQSNHANIFYSVLTIKMQNKDKGLYTCRVRSGPSFKSVNTS VHIYDKAFITVKHRKQVLETVAGKRSYRLSMKVKAFFSPEVWVWKDGLPATEKSARYLTRGYSIIKDVTEED AGNYTILLSIKQSNVFNLTATLIVNVKPIYKAVSSFPDPALYPLGSRQILTCTAYGIPQPTIKWFWHPCNHNH SEARCDFCSNNEESFILDADSNMGNRIESITQRMAIEGKNKMASTLVVADSRISGIYICIASNKVGTVGRNISFYI TDVPNGFHVNLKMPTEGEDLKLCTVKNFLYRDVTWILLRTVNNRTMHYSISKQKMAITKEHSITLNLTIMNVS LQDSGTACRARNVYTGEEILQKKEITIRDQEAPYLLRNLSHDHTVAISSSTLDCHANGVPEPQITWFKNNHKIQ QEPGILGPGSSTLFIERVTEEDGVYHCKATNQGKSVESAYLTVQGTSDKSNVDDIEGRMDEPKSCDKTHT CPPCPAPPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYN STYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVK GFYPSDIAVEWESNGQPENNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNVFCSSVMHEALHNHYTQKSL SLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only