

Recombinant Rat VEGF-A, VEGF164

Catalog No: #AP76130



Package Size: #AP76130-1 10ug #AP76130-2 50ug #AP76130-3 500ug #AP76130-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Rat VEGF-A, VEGF164
Host Species	Yeast
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Rat
Accession No.	P16612-2
Uniprot	P16612
GeneID	83785;
Calculated MW	19.2kDa
Target Sequence	APTTEGEQKTHEVVKFMDVYQRSYCRPIETLVDIFQEYPDEIEYIFKPSCVPLMRCAGCCNDEALECVPTSESN VTMQIMRIKPHQSQHIGEMSFLQHSRCECRPKKDRTPENHCEPCSERRKHLFVQDPQTCKCCKNTDSRC KARQLELNERTCRCDKPRR
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only