

Recombinant Vibrio cholerae Toxin B

Catalog No: #AP76134



Package Size: #AP76134-1 10ug #AP76134-2 50ug #AP76134-3 500ug #AP76134-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Vibrio cholerae Toxin B
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Vibrio cholerae
Accession No.	P01556
Uniprot	P01556
GeneID	2613962;
Calculated MW	11.8kDa
Target Sequence	MTPQNITDLCAEYHNTQIYTLNDKIFSYPESLAGKREMAITFKNGAIFQVEVPGSQHIDSQKKAIERMKDTRLRIAYL TEAKVEKLCVWNNKTPHAIAAISMAN
Formulation	Lyophilized from a 0.2 µm filtered solution of 50mM Tris, 200uM NaCl, pH8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only