

# Recombinant Human Sulfotransferase,SULT1C4,SULT1C2

Catalog No: #AP76136

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP76136-1 10ug #AP76136-2 50ug #AP76136-3 500ug #AP76136-4 1mg

## Description

|                 |  |
|-----------------|--|
| Product Name    | Recombinant Human Sulfotransferase,SULT1C4,SULT1C2   |
| Host Species    | E.coli   |
| Purification    | Greater than 95% as determined by reducing SDS-PAGE.   |
| Specificity     | Human  |
| Accession No.   | Q6PD90   |
| Uniprot         | Q6PD90   |
| Calculated MW   | 14.1kDa  |
| Target Sequence | MGSSHHHHHSSGLVPRGSHMMALHDMEDFTFDGTRKLSVNYVKILQPTDTCDIWDKIWNFQAKPDDLIS<br>TYPKAGTTWTQEIVELIQNEGDVEKSKRAPTHQRFPFLEMKIPSLGSGEYK   |
| Formulation     | Lyophilized from a 0.2 µm filtered solution of 4mM HCl.  |
| Storage         | Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months. |

Note: This product is for in vitro research use only