

# Recombinant Human $\beta$ -1,3-Glucuronyltransferase 3,B3GAT3 (C-6His)

Catalog No: #AP76138

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP76138-1 10ug #AP76138-2 50ug #AP76138-3 500ug #AP76138-4 1mg

## Description

Product Name	Recombinant Human $\beta$ -1,3-Glucuronyltransferase 3,B3GAT3 (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O94766
Uniprot	O94766
GeneID	26229;
Calculated MW	30.4kDa
Target Sequence	MEALPTIYVVTPTYARLVQKAEVRLSQTLSLVPRLHWLLVEDAEGPTPLVSGLLAASGLLFTHLVLTTPKAQRL REGEPGWVHPRGVQRNKALDWLRGRGGAVGGEKDPPTGQGVVYFADDNTYSRELFEEMRWTRGVS VWPVGLVGGLRFEGPQVQDGRVVGFTAWEPSRPFVDMAGFAVALPLLLDKPNAQFDSTAPRGHLESSLL SHLVDPKDLEPRAANCTRLVWHTRTEKPKMKQEEQLQRQGRGSDPAIEVLEHHHHHH
Formulation	Supplied as a 0.2 $\mu$ m filtered solution of 20mM Tris,150mM NaC,2mM EDTA,20% Glycerol,pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only