

# Recombinant Human $\delta(3,5)$ - $\delta(2,4)$ -Dienoyl-CoA Isomerase, Mitochondrial, ECH1



Catalog No: #AP76140

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP76140-1 10ug #AP76140-2 50ug #AP76140-3 500ug #AP76140-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human $\delta(3,5)$ - $\delta(2,4)$ -Dienoyl-CoA Isomerase, Mitochondrial, ECH1
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q13011
Uniprot	Q13011
GeneID	1891;
Calculated MW	34.45kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMTGSSAQEEASGVALGEAPDHSYESLRVTSAQKHVLHVQLNRPNKRNAMN KVFWREMVECFNKISRDCRAVVISGAGKMFTAGIDLMDMASDILQPKGDDVARISWYLRDIITRYQETFNVI ERCPKPVIAAVHGGCIGGGVDLVTACDIRYCAQDAFFQVKEVDVGLAADVGTQLRQLPKVIGNQSLVNELAF RKMMADEALGSLVSRVFPDKEVMLDAALALAAEISSKSPVAVQSTKVNLLYSRDHLSVAESLNYVASWNMSM LQTQDLVKSVQATTENKELKTVTFSKL
Formulation	Supplied as a 0.2 $\mu$ m filtered solution of 20mM Tris, 100mM NaCl, 10% Glycerol, pH 8.0.
Storage	Store at $-20^{\circ}\text{C}$ , stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only